

DEPARTMENT OF ELECTRONICS AND TELECOMMUNICATION ENGINEERING

BM 3121 MEDICAL IMAGING

Ultrasound Quality Assurance

Name:

U.G.C. Jayasankha

Index Number: 170259P

Contents

| 1 | | roduction to Diagnostic Ultrasound | 2 |
|---|-----|--|----|
| | 1.1 | | 2 |
| | 1.2 | Importance of Diagnostic Ultrasound | |
| | 1.3 | Components and Implementation of Diagnostic Ultrasound | - |
| 2 | Qua | ality Assurance of Diagnostic Ultrasound | 6 |
| 3 | Pha | antoms for Diagnostic Ultrasound QA | 8 |
| | 3.1 | TMM based Phantoms | 8 |
| | 3.2 | Water based phantoms | Ć |
| 4 | Ma | in Quality Assurance Methods of Diagnostic Ultrasound | 10 |
| | 4.1 | Visual Checklist | 10 |
| | | 4.1.1 Infection control and scanner damage | 10 |
| | | 4.1.2 Basic scanner and transducer testing | 10 |
| | | 4.1.3 Further scanner and transducer testing | 10 |
| | 4.2 | System Sensitivity test | 12 |
| | | 4.2.1 Maximum Depth of Penetration/Visualization | 12 |
| | | 4.2.2 Speckles, artefacts and Noise | 13 |
| | 4.3 | Uniformity test | 14 |
| | | 4.3.1 Horizontal Uniformity | 14 |
| | | 4.3.2 Vertical Uniformity | 14 |
| | 4.4 | Spatial Accuracy test | 15 |
| | | 4.4.1 Horizontal and vertical Distance Accuracy | 15 |
| | | 4.4.2 Volume Measurement | 16 |
| | 4.5 | Spatial resolution test | 16 |
| | 4.6 | Contrast Resolution test | 17 |
| | | 4.6.1 High Contrast Detectability test - Anechoic objects | 17 |
| | | 4.6.2 Low Contrast Object Detectability test - Soft tissue | 17 |
| 5 | Rec | cent advancements in Ultrasound Quality Assurance | 19 |
| 6 | Ref | ference | 19 |

1 Introduction to Diagnostic Ultrasound

1.1 Basics of Diagnostic Ultrasound

Medical Ultrasound or Diagnostic-Sonography is a medical procedure of using high frequency sound waves to image inside of the body. Ultrasound imaging was invented in 1956. It can be used to display dynamic visual images in most parts of the body such as tendons, muscles, joints, blood vessels, and internal organs. Ultrasound imaging is the most popular imaging modality after projection x-ray imaging. Ultrasound imaging is predominantly used to image and diagnose pregnant mothers, which is known as "obstetric ultrasound" since ultrasound does not use ionizing radiations for imaging, thus no harm to the fetus. Ultrasound is also highly used to image soft tissues such as breast, tendons, thyroid, parathyroid glands and testes.



Figure 1: Ultrasound Imaging of a pregnant woman

Physics of the ultrasound imaging is the echoing nature of sound waves. When the sound waves travel through a medium it attenuates. When the sound waves hit a boundary part of the wave transmits to the second medium and other part reflects back to the first medium. When this reflecting wave is detected it is proof for existing boundary. And using the time difference the distance to the boundary can be calculated if the speed of sound wave is known. Several reflection wave fronts will arrive at the receiver if there are several boundaries.

There are 4 modes in ultrasound imaging.

- i. **A-mode / Amplitude mode** the simplest type of ultrasound. A single transducer scans a line through the body with the echoes plotted on screen as a function of depth.
- ii. **B-mode** / **Brightness mode** is the display of 2D map of B-Mode data, and is the most common form of ultrasound imaging. Unlike A-Mode, B-Mode is based on brightness with the absence of vertical spikes. The brightness depends upon the amplitude or intensity of the echo.
- iii. M-mode / Motion mode a rapid sequence of B-mode scans whose images follow each other in sequence on screen enables doctors to see and measure range of motion, as the organ boundaries that produce reflections move relative to the probe.
- iv. **Doppler mode** This mode makes use of the Doppler effect in measuring and visualizing blood flow. Doppler sonography play important role in medicine. Sonography can be enhanced with Doppler measurements, which employ the Doppler effect to assess whether structures (usually blood) are moving towards or away from the probe, and its relative velocity.




Figure 2: Modes of Ultrasound

1.2 Importance of Diagnostic Ultrasound

Ultrasound scanning is ranked second after x-ray imaging as most used imaging modality. There are several reasons for that.

- 1. Ionizing radiation is not used
 - X-ray, CT, PET modalities of imaging use ionizing radiation which is harmful for the organs. But in Ultrasound imaging; ionizing radiation is not used. Therefore it is not harmful for the organs.
- 2. Inexpensive equipment & low cost scans
 When compared with other imaging modalities, Ultrasound scanners are less
 complicated. Therefore the system is inexpensive and maintenance cost is low.
 People can afford Ultrasound imaging modality easily.
- 3. High temporal resolution

 Ultrasound can be operated in hi

Ultrasound can be operated in higher frame rates. Thus images can be acquired in short periods of time. Movements of organs such as moving baby, beating heart can be seen because of high frame rate.

4. Flow velocity measurement availability
Using Doppler effect, velocity of moving fluids in the body can be calculated.
This can be useful to identify flow direction. Fluid moving patterns can be modeled.

5. Portability

Ultrasound scanner are light and smaller when compared with other modalities. It does not need an dedicated imaging room. It can be used in any ward and can be transported somewhere else easily.

But there are several disadvantages also.

- 1. Ultrasound imaging is not usable in bony or airy regions of body. Ex: Lungs
- 2. Field of view is limited. Can not get clear, accurate images as MRI or X-ray
- 3. Quality of imaging depends on the operator.
- 4. Ultrasound images are hard to interpret. Experience is a must to interpret ultrasound images.
- 5. Artefacts are common.

1.3 Components and Implementation of Diagnostic Ultrasound



Figure 3: Main components of Diagnostic Ultrasound

All the main components of US system are shown in above figure. After applying gel on the imaging area of the body, using the transducer imaging is done. This is less complicated than other modalities.




Figure 4: Imaging system block diagram

2 Quality Assurance of Diagnostic Ultrasound

Ultrasound is a widely used imaging modality when compared with other imaging modalities since they are cheaper and harmless. However quality assurance of Diagnostic ultrasound is not a popular topic as quality assurance of other imaging modalities. Other imaging modalities are very complex and expensive equipment. Thus quality assurance methods are followed strictly. But due to some reasons people tend to neglect quality assurance methods of ultrasound. Here are the reasons for that given by people who refuses QA measures,

- i. Ultrasound is considered as a safe imaging modality. Thus QA measures are not needed
- ii. Ultrasound Imaging systems are stable. Therefore QA measures are not necessary.
- iii. There are no dedicated operator or engineer for Ultrasound systems. Even the systems are being moved to different sections. Thus nobody knows how, what and when to test.
- iv. No well established rules and regulations.



Figure 5: Movable Ultrasound scanners

But these reasons are not accurate. Even though ultrasound scanners were used more than 50 years as a safe modality it had been revolutionary developed through past few years. In the beginning it was a simple system but now it is a compact system with various components. Due to competitive nature between vendors of ultrasound systems they try to overcome other using new technologies. Therefore most of the new ultrasound scanners use high power output. Even though it is not proved, there are evidence that high power outputs can harm the subject. Therefor quality assurance methods and rules and regulations are needed to limit exposure time, minimize scan times, reduce repeat rates and assure that scan protocols employ proper scanner settings. Even though Ultrasound systems are stable than other modalities, faults can occur gradually. A quality assurance measures can identify these faults or damages at the beginning. So it can be repaired with a low cost.

Not every Ultrasound scanner is not expected to exhibit failures but ultrasound equipment is not immune to failures. Therefore Quality assurance measures are expected to be followed. Otherwise failures may go unnoticed and cause sever damages to the system and diagnosis may be wrong.

As in the other imaging modalities **phantoms** are used for quality assurance tests of Diagnostic Ultrasound.

Following are the Quality assurance test and methods for Diagnostic Ultrasound

- 1. System Sensitivity test
 - a. Maximum Depth of Penetration
 - b. Speckles, artefacts and Noise
- 2. Spatial Accuracy test
 - a. Volume Measurement
 - b. Horizontal and vertical Distance Accuracy
- 3. Uniformity test
 - a. Horizontal Uniformity
 - b. Vertical Uniformity
- 4. Contrast Resolution test
 - a. High Contrast Detctability test
 - b. Low Contrast Object Detctability test

- 5. Spatial resolution test
 - a. Horizontal resolution
 - b. Vertical Resolution

3 Phantoms for Diagnostic Ultrasound QA

Phantoms are used for quality assurance purposes to measure relevant parameters. These phantoms are specially designed for Ultrasound to represent the human body and mimic the human nature. Following important features should be similar in phantom and the human body.

- i. Attenuation Coefficient
- ii. Speed of Sound
- iii. Back-scatter Coeffucient
- iv. Mechanical Properties
- v. Thermal Properties
- vi. Elasticity

3.1 TMM based Phantoms

There are special substances which have above features similar to human tissues. Those substances are known as Tissue Mimicking Materials(TMM). They are used to manufacture Ultrasound phantoms.

- i. Urethane
- ii. Agar
- iii. Zerdine
- iv. TPE

These are solid non-water based materials. Advantages of these TMM are,

- Not subject to desiccation
- No need for scanning window
- Produce tissue-like backscatter

Disadvantages are,

- Velocity is lower than tissue
- Attenuation is not proportional to frequency.
- Surface easily gets damaged if not cleaned




Figure 6: TMM based phantom

3.2 Water based phantoms

Water based phantoms have some advantages over TMM based phantoms,

- Velocity is similar to tissue
- Attenuation is proportional to frequency.
- Back-scattering property

But there are few disadvantages also,

- Subject to desiccation
- Should be kept in containers
- Scanning window is required



Figure 7: Water based phantom

4 Main Quality Assurance Methods of Diagnostic Ultrasound

This section describes the main quality assurance methods used in diagnostic ultrasound.

4.1 Visual Checklist

This is the basic method to find flaws in the system and it should be done in every week. Following list of equipment is inspected visually to find any missing part or any malfunction. This verifies whether the US system working properly electrically and mechanically.

Visual checklist is divided into three sections

4.1.1 Infection control and scanner damage

4.1.2 Basic scanner and transducer testing

4.1.3 Further scanner and transducer testing

If any of this fails, it should be repaired or replaced as soon as possible even before carrying out other quality assurance protocols.

| Identifier | Task | Frequency |
|------------|--|-------------|
| 1.1 | Clean ultrasound gel and body fluids from scanner console, transducers and cables after every patient | Multi-daily |
| 1.2 | Ensure transducers are stored securely when not in use | Multi-daily |
| 1.3 | Ensure transducer cables are properly stowed and not at risk of being run over by the scanner wheels | Multi-daily |
| 1.4 | Clean monitors of dust, gel, etc. | Daily |
| 1.5 | Check the correct operation of the main scanner controls | Daily |
| 1.6 | Inspect the transducers used during the session for damage | Daily |
| 1.7 | Inspect switches, knobs and other controls for damage | Weekly |
| 1.8 | Inspect probe cables, the mains cable and plug and other cables, e.g. network and printer cables, for damage | Weekly |
| 1.9 | Inspect the ultrasound system for damage such as cracks and dents | Weekly |
| 1.10 | Test of brake and wheel function | Weekly |
| 1.11 | Check air filters for dust and fluff | Weekly |

Figure 8: Visual checklist step 1

| Identifier | Task | Frequency |
|------------|---|-----------|
| 2.1 | Check that the video monitor's brightness and contrast controls have been appropriately adjusted | Daily |
| 2.2 | Check that the greyscale bar is fully displayed | Daily |
| 2.3 | Check grey levels on the ultrasound system monitor is matched by grey levels on PACS system monitors | Daily |
| 2.4 | Inspect reverberation images for shadows and streaks caused by transducer dropout for all transducers in use that day | Daily |

Figure 9: Visual checklist step 2

| Identifier | Task | Frequency |
|------------|--|--|
| 3.1 | Air reverberation pattern (sensitivity) | Monthly |
| 3.2 | Element dropout test | If fault in 2.4. Monthly for phased arrays |
| 3.3 | Electronic noise assessment (B-mode, Doppler and Colour Flow) | Monthly |
| 3.4 | Date of electrical safety testing | Monthly |

Figure 10: Visual checklist step 3

4.2 System Sensitivity test

4.2.1 Maximum Depth of Penetration/Visualization

Maximum visual depth or maximum penetration depth is the maximum depth which can be seen through image. This is a axial measure of sensitivity. It is considered by as a overall check of the integrity of the US system.



Figure 11: Maximum Depth of Penetration

4.2.2 Speckles, artefacts and Noise

Speckles - Abnormalities arose from the destructive and constructive nature of sound waves. This is due to the reflection of sound waves hitting small bodies.



Figure 12: Speckles in Diagnostic Ultrasound

Artefacts - Abnormalities that does not have any relation with the actual bodies/structures. Speckles are also artefacts in general. **Noise** - Noise is the unwanted



Figure 13: Artifacts in Diagnostic Ultrasound

additional parts of the signal/image. Imperfections of the transducer or any other part of the system can generate noises. To achieve a great quality image, this noise part should be suppressed. **SNR Measurement** is the parameter to identify the noise level. SNR ratio decreases with the depth increases.




Figure 14: SNR vs Depth

4.3 Uniformity test

Uniformity test is considered as the most important and useful test. In an ideal image,

- No evidence of elements dropout
- No vertical shadows
- No loss of sensitivity near edges

4.3.1 Horizontal Uniformity

Horizontal non-uniformity is caused by malfunctioning transducer unit/element in the array. This can be identified by using a phantom. Following figure shows a such situation.




Figure 15: Horizontal Uniformity defects

4.3.2 Vertical Uniformity

Vertical non-uniformity occurs due to the attenuation of the sound waves. However, Time Gain Compensation (TGC) is used to produce a vertically uniform output. When Time Gain Compensation has any fault vertical non-uniformity can occur.




Figure 16: Vertical Uniformity defects

4.4 Spatial Accuracy test

4.4.1 Horizontal and vertical Distance Accuracy

Using a known phantom which has predefined structure ultrasound image is obtained. Obtained image is compared with the original image and the spatial accuracy of the system can be evaluated. If they differ a lot more than accepted values given by the phantom vendor, system should undergo a repair.




Figure 17: Horizontal and vertical Distance Accuracy test

4.4.2 Volume Measurement

In 3D ultrasound, volumetric information heavily depends on accuracy of the linear measurements in all three dimensions. Therefore it is necessary to evaluate this using the same procedure described above using 3D targets.

4.5 Spatial resolution test

Minimum spatial resolution is the smallest separation that can be seen on the image. Evaluation of spatial resolution might also be performed by using pixel data across a bright target, then plotting to determine FWHM. Both are affected by the gain and gray scale processing settings Following diagram shows a method of spatial resolution measurement.




Figure 18: Spatial resolution test

Spatial resolution is,

- Frequency/transducer dependent
- Depth dependent
- specified at a frequency and a depth of phantom

4.6 Contrast Resolution test

4.6.1 High Contrast Detectability test - Anechoic objects

Ability to detect anechoic (weakly echogenic objects which reflects sound waves weakly) objects in off-axis is known as contrast resolution. Example for an anechoic object is a cyst. Weakly echogenic objects such as blood vessels also get affected by this. Clutter can be caused by several sources such as S side lobes, quantization clutter and grating lobes.

Contrast Resolution in this case is defined in terms of Clutter Energy to Total Energy Ratio (CTR). This is measured by taking the average pixel intensity of a cyst phantom image compared to the surrounding.




Figure 19: High Contrast Detectability test

4.6.2 Low Contrast Object Detectability test - Soft tissue

Low contrast resolution/ soft tissue resolution refers to the ability to distinguish echogenicity differences between neighboring soft tissue regions. As example,

- a hyper-echoic or hypo-echoic lesion in normal soft tissue (liver parenchyma)
- different types of tissue side by side(liver/bowel/kidney/liver)

Speckles in ultrasound images makes detecting subtle echogenicity variations difficult. Tissue contrast resolution can be increased by reducing speckle variance or by reducing speckle size.




Figure 20: Low Contrast Detectability test

5 Recent advancements in Ultrasound Quality Assurance

Diagnostic Ultrasound systems are simple when compared with other imaging modalities. They have developed very well and can be assumed as its peak now. However, currently many research efforts are carried out to discover new ways to enhance the level of quality assurance of ultrasound imaging. Some of the quality assurance researches related to Ultrasound imaging quality assurance can be listed as follows.

- Image quality improvements
- Advancements in volumetric ultrasound
- Ultrasound in every physician's pocket
- Sonoelastography

6 Reference

- Dudley, N., Russell, S., Ward, B., Hoskins, P., & BMUS QA Working Party. (2014, February). BMUS guidelines for the regular quality assurance testing of ultrasound scanners by sonographers. Retrieved May 17, 2020, from https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4760519/
- 2. "Medical ultrasound," Wikipedia, 17-May-2020. [Online]. Available: https://en.wikipedia.org/[Accessed: 18-May-2020].
- 3. A. R. Webb, Introduction to biomedical imaging. Hoboken, N. J.: John Wiley & Sons, 2006.
- 4. Evan J. Boote, Ph.D. University of Missouri-Columbia "Current Ultrasound Quality Control Recommendations and Techniques"
- 5. Kutay F. Üstüner and Gregory L. Holley Siemens Medical Solutions USA, Inc. Ultrasound Division Mountain View, "CA Ultrasound Imaging System Performance Assessment"
- 6. "Home," RSS. [Online]. Available: http://www.aapm.org/. [Accessed: 18-May-2020].

Appendix

List of Figures

| 1 | Ultrasound Imaging of a pregnant woman |
|----|--|
| 2 | Modes of Ultrasound |
| 3 | Main components of Diagnostic Ultrasound |
| 4 | Imaging system block diagram |
| 5 | Movable Ultrasound scanners |
| 6 | TMM based phantom |
| 7 | Water based phantom |
| 8 | Visual checklist step 1 |
| 9 | Visual checklist step 2 |
| 10 | Visual checklist step 3 |
| 11 | Maximum Depth of Penetration |
| 12 | Speckles in Diagnostic Ultrasound |
| 13 | Artifacts in Diagnostic Ultrasound |
| 14 | SNR vs Depth |
| 15 | Horizontal Uniformity defects |
| 16 | Vertical Uniformity defects |
| 17 | Horizontal and vertical Distance Accuracy test |
| 18 | Spatial resolution test |
| 19 | High Contrast Detectability test |
| 20 | Low Contrast Detectability test |